


COVID-19 TESTING

Nasal Swab Self-Collection Instructions


A sample has been requested by your employer for testing. This test will help identify whether you have a current COVID-19 infection. Please follow the instructions below carefully to safely self-collect and ship your test sample to the Biodesix laboratory.

SAMPLE COLLECTION AND SHIPPING INSTRUCTIONS

- 1. Label your sample tube.** Using the provided labels, place your first and last name, date of birth, and date of collection on the label, then stick the label on the sample tube.




- 2. Open nasal swab.** Remove from sterile packaging by pulling the two ends of the wrapper apart (like a bandage).




NOTE: Be careful not to touch the soft tip

- 3. Swab nose to collect sample.** Gently place the entire soft tip of the swab into one nostril until you feel resistance. Rub it in a circle 4 times inside your nostril while placing gentle pressure with your finger on the soft tip from the outside of the nostril to obtain the best sample for testing. Repeat by placing the same swab into the other nostril and circle 4 more times.



IMPORTANT: Both nostrils must be sampled to get a proper specimen.

- 4. Place the swab into the tube.** Carefully lower the swab, tip first, into the provided tube. Once the tip is at the bottom, break the swab handle at the top of the tube by bending it back and forth. Then, screw the cap on tightly.



CAUTION: If any of the liquid in the sample tube is spilled, do NOT use bleach to clean up the liquid. Clean the spill with alcohol, followed by water.

- 5. Package and ship** according to the instructions below.

PACKAGING AND SHIPPING INSTRUCTIONS

1. Confirm that the label is applied to the sample tube. Samples received at the laboratory without a label will need to be re-collected.



2. Complete the Sample Information card and place it in the return shipment box.

PLEASE RETURN THIS CARD WITH YOUR SAMPLE

biodesix


COVID-19 TESTING
Sample Information
 Please complete the information below and return with your sample to help the laboratory contact your sample to your testing program. If the laboratory has specific questions about your sample, we will use the contact information below to reach out to you.

Employer Name: Employer, Inc
 First Name: John
 Last Name: Smith
 Date of Birth (MM/DD/YYYY): 10/23/1967
 Contact email: John.smith@employer.com
 Contact phone number: ###-###-####


For questions, please contact covidtesting@biodesix.com or call 1.833.578.8559.

Employer Name: Employer, Inc
 First Name: John
 Last Name: Smith
 Date of Birth (MM/DD/YYYY): 10/23/1967
 Contact email: John.smith@employer.com
 Contact phone number: ###-###-####

3. Place sample tube into the biohazard bag and place into the return shipment box with the sample information card.



4. Place return shipment box into the pre-labeled FedEx UN 3373 shipping envelope.




5. Contact FedEx immediately at 1.800.463.3339 to locate the nearest drop box or schedule a pickup for overnight shipment.



Biodesix, Inc.
 2970 Wilderness Place Suite 100
 Boulder, CO 80301

Made in the USA
www.biodesix.com



Do not use if the product sterile barrier system or its packaging is compromised.